

CORRIGENDUM01


Tender Ref. No.:- BHU/ZOO/DS/EQUIP/2020-21/01 **Dated:** 24-11-2020

Tender Id :- 2020_BHU_598800_1

Tender Titled :- Supply of Chemi-Gel Doc Unit & Refrigerated Centrifuge

Corrigendum Titled : - Amended specifications of Refrigerated centrifuge

Old specifications	New specifications
<ol style="list-style-type: none">1. Swing out rotor of 4x 750 ml to universally fit tube/bottle/Plates using universal adapter with maximum of 4,500 rpm or more and 5,000 x g or more<ol style="list-style-type: none">a. Adapters for Tubes capacity: 45 x 15 mL or more, 16 x 50 mL or more, 4 x 250 mL bottle and 4 x MTP/DWP2. Fixed angle rotor for 6 x 50 mL with aerosol tight lid for 50 mL conical tubes/18 ml oakridge tube (round bottom tubes) with maximum of 14,000rpm and 21,000 x g or more<ol style="list-style-type: none">a. Adapters for 15 mL tubes must be providedb. Adaptors for 18 mL oakridge tube (round bottom tubes) must be provided3. Fixed angle rotor 30 x 1.5 / 2 mL with aerosol tight lid with maximum of 13,000 rpm and 20,000 x g or more.<ol style="list-style-type: none">a. Adaptors for 0.5/0.6 mL tubes to be provided4. A fixed angle rotor, 24 x 1.5 / 2 mL with aerosol tight lid with maximum of 12,000 rpm and 15,000 x g or more, containing module with a flexible installation base for quick running of spin columns at room temperature as per the demands of lab protocols.	<ol style="list-style-type: none">1. Swing out rotor of 4x 750 ml to universally fit tube/bottle/Plates using universal adapter with maximum of 4,500 rpm or more and 5,000 x g or more<ol style="list-style-type: none">a. Adapters for Tubes capacity: 45 x 15 mL or more, 16 x 50 mL or more, 4 x 750 mL (if required), 4 x 250 mL bottle and 4 x MTP/DWP2. Fixed angle rotor for 6 x 50 mL with aerosol tight lid for 50 mL conical tubes and ~50 ml oakridge tube (round bottom tubes) with maximum of 14,000rpm and 21,000 x g or more<ol style="list-style-type: none">a. Adapters for 15 mL tubes must be providedb. Adaptors for ~50 mL oakridge tube (round bottom tubes) must be provided3. Fixed angle rotor 30 x 1.5 / 2 mL with aerosol tight lid with maximum of 13,000 rpm and 20,000 x g or more.<ol style="list-style-type: none">a. Adaptors for 0.5/0.6 mL tubes to be provided4. A fixed angle rotor, 24 x 1.5 / 2 mL with aerosol tight lid with maximum of 12,000 rpm and 15,000 x g or more, containing module with acceleration/ deceleration of maximum 15 sec with continuous running capability of at least 15 sec to 9h, having a flexible installation base for quick running of spin columns at room temperature as per the demands of lab protocols. <p>Optional:</p> <ol style="list-style-type: none">1. An 18-place aerosol tight rotor containing device designed for 1.5/2 mL tubes with speeds of at least 16,000 x g in a refrigerated set up, having maximum acceleration/deceleration time of 15 sec and continuous run function of at least 15 sec to 9h.


(Devanjan Sinha)
Dr. Devanjan Sinha
Assistant Professor
Department of Zoology
Institute of Science
Banaras Hindu University
Varanasi-221005, India